

Shimadzu LCMSsolution

for
LCMS-2010 / LCMS-QP8000 α

Operation Guide

Read the instruction manual thoroughly before you use the product. Keep this instruction manual for future reference.



Shimadzu Corporation

Analytical & Measuring Instruments Division

Kyoto, Japan

Copyright Shimadzu Corporation 2003. All rights are reserved, including those to reproduce this publication or parts thereof in any form, without the express written permission of Shimadzu Corporation.

Information in this publication is subject to change without notice and does not represent a commitment on the part of the vendor.

Any errors or omissions which may have occurred in this publication will be corrected as soon as possible, but not necessarily immediately upon detection.

Shimadzu Corporation makes no claims or warranty concerning the influence of the application of this publication on the operation of the machine.

Note that Shimadzu does not have any obligation concerning the effects resulting from the application of the contents of this manual.

Shimadzu Corporation will continue to supply original Shimadzu maintenance parts associated with a given product for a period up to 7 years from the cessation of production of that product.

Please be informed that after this period of time, Shimadzu Corporation cannot guarantee supply of such original maintenance parts. However, following discontinuation of a product, the period of availability of maintenance parts which have been produced on a subcontract basis is up to the discretion of the concerned sub-contracting company.

Microsoft® Windows®, Windows NT® and MS-Excel® are registered trademarks of Microsoft Corporation. Microsoft® Windows® 2000 and MS-Word are products of Microsoft Corporation.

Adobe®, Acrobat® and Distiller® are registered trademarks of Adobe Systems Incorporated. Acrobat Reader™ and Acrobat PDFWriter™ are products of Adobe Systems Incorporated.

Introduction

Thank you very much for purchasing the LCMSsolution software for Shimadzu liquid chromatography / mass spectrometry workstations (hereafter called "LCMSsolution").

LCMSsolution allows you to control the liquid chromatograph (hereafter called "LC") and the Mass Spectrometer (hereafter called "MS") from your personal computer, acquire chromatograms and other different kinds of data, and reanalyze the acquired data under different parameters on your personal computer.

This manual is the tutorial in the most simplified analysis procedure using LCMSsolution which helps you to catch more knowledge in other volumes or further actual operations.

The "Operation manual" and "Administration manual" are attached as separate volumes.

The Operation manual has been put together in order to familiarize you with the basic knowledge required to operate LCMSsolution. Be sure to read it thoroughly before using this software. After reading the manual, keep it in a safe place so that it can be accessed whenever necessary.

The Administration manual covers the information useful for system administration such as the support features for GLP/GMP or FDA 21CFR Part11, a set of regulations for electronic records and electronic signature. For more information on the functions of LCMSsolution, refer to this on-line manual.

This manual assumes that the reader is knowledgeable of basic operations of Windows®2000. For the operation of Windows®2000, refer to the instruction manual that comes with that product.

This manual sometimes explains commonly for LabSolutions series. And some explanations may use the drawings come from sister products like LCsolution, if it does not cause misunderstanding in the range of explanations.

Using the instruction manual




Kinds of instruction manuals


The LCMSsolution package contains the following information that describes the operational procedures and functions.

Name	Media	Description
Operation guide for LCMSsolution	Printed Document	Provides tutorial on mostly basic analysis procedure using LCMSsolution.
Operation manual for LCMSsolution	Printed Document	Explains the operational procedures for data acquisition and analysis using LCMSsolution.
Administration manual for LCMSsolution	Printed Document	Explains the operational procedures and basic idea of system administration and data management using LCMSsolution.
On-line help	LCMSsolution program	Provides detailed information on parameters and setting ranges. This is accessible from the Help menu in LCMSsolution. (For using the on-line help, refer to section "14.1.1 Using Help" in the Operation manual.)
Operation guide for LCMSsolution (PDF version)	CD-ROM disk for installation	Provides the operation guide volume of the instruction manual as a PDF file so that it can be viewed on your personal computer. The general table of contents is available, including other instruction manuals (PDF versions). It allows you to use each instruction manual via the hyperlink.
Operation manual for LCMSsolution (PDF version)	CD-ROM disk for installation	Provides the operation volume of the instruction manual as a PDF file so that it can be viewed on your personal computer. It is accessible from the Help menu in LCMSsolution. (For using this PDF, refer to section "14.1.2 Using the Online Manual" in the Operation manual.)
Administration manual for LCMSsolution (PDF version)	CD-ROM disk for installation	Provides the administration volume of the instruction manual as a PDF file so that it can be referred to on-line whenever operations related to system administration are needed. The general table of contents is available, including all the instruction manuals (PDF versions). It allows you to use each instruction manual via the hyperlink.

Legends for instruction manual

This manual uses the following legends:

Legend	Meaning
	Shows additional informations around the topic.
	Points the reference informations.
	Gives you tips.
< >	Shows a window or view name; e.g., <Data Acquisition> window or <Method> view.
[]	Shows a parameter, tab, column, cell, bar name, menu command, that can be selected from the menu bar.
[]-[] command	Shows a sequence of selecting the menu in the first [] and then selecting the command in the second []. For example, [File]-[Print] command means that you should click on the File menu and then select the Print command from the displayed list of commands.



Contents

1 Making Preparations for Analysis	1
1.1 Basics of LCMSsolution.....	1
1.2 Starting the LCMSsolution.....	3
2 Qualitative Processing (Single-run Analysis)	6
2.1 Creating a new method file.....	6
2.2 Setting the LC parameters.....	7
2.2.1 Detecting the auto sampler rack	7
2.2.2 Setting the LC parameters	7
2.3 Setting the MS parameters.....	9
2.4 Starting the operation of the instrument	11
2.4.1 Starting the control of the MS unit	11
2.4.2 Starting the operation of the LC unit	12
2.4.3 Selecting a graph to be displayed in the <Chromatogram> view	13
2.5 Acquiring data through a single-run analysis.....	14
2.6 Performing qualitative processing on <MS Data Analysis>	15
2.6.1 Starting the <MS Data Analysis>	15
2.6.2 Displaying a mass spectrum	16
2.6.3 Displaying a mass chromatogram	18
2.7 Performing peak integration (peak detection).....	19
2.8 Printing out the analysis result.....	21
2.8.1 Printing out a “Graph Image”	21
2.8.2 Selecting a layout for printing	23
3 Quantitative Processing (Batch Analysis)	25
3.1 Creating a “Compound Table”	25
3.1.1 Setting the quantitative parameters in <MS Data Analysis>	25
3.1.2 Creating a “Compound Table”	26
3.2 Creating a SIM Table.....	28
3.3 Creating a “Batch Table”	29
3.4 Making a batch analysis	32


4 Data Analysis	36
4.1 Checking a “Calibration Curve”	36
4.2 Checking the quantitative calculation result of an unknown sample.....	38
4.3 Loading a batch file to the “Quant Browser”	39
4.3.1 Displaying the quantitative result from the batch file	39
4.3.2 Setting the integration parameters again to retry peak integration	40
4.4 Printing out a summary report from the Quant Browser	43
4.5 Using the Data Browser.....	44
5 Exiting the LCMSsolution	46
5.1 Existing the LCMSsolution.....	46
Index	





1

Making Preparations for Analysis

1.1 Basics of LCMSsolution

■ <LCMSsolution Launcher> - [Operation] menu icon



No.	Icon	Name	Description
1		Analysis	Starts the application for configuring and controlling the system and making a single-run or batch analysis. (Starts <LCMS Analysis> in the Online mode)
2		Offline Editor	Starts the application for editing any method file or batch file not in use during the analysis. (Starts <LCMS Analysis> in the Offline mode)
3		Postrun	Starts the application for loading the acquired analysis data to create a calibration curve or perform data processing.
4		Browser	Starts the application for browsing multiple analysis data together or analyzing data together.


Files used in LCMSsolution

Extension	Name	Description
.lcm	Method File	Analysis condition, Data processing conditions, QA/QC settings, calibration curve information, and system configuration
.lcr	Report Format File	Report formats
.lcb	Batch File	Batch tables and batch settings
.lcd	Data File	Chromatograms, mass spectrums, peak tables, identification/quantitation results, report format, tuning results, methods, and batch table

 [Admin Manual]: [“4.1 Important File Concepts for Operation”](#)

Data structure in LCMSsolution

The data in the LCMSsolution is retained in data files, consisting various types of records and parameters such as the system configuration, fine-tuning result, system conditions, and analysis conditions that have been used to acquire and analyze data. This structure enables you to browse each data file for monitoring conditions and analysis parameters, thereby ensuring the traceability of data. This means that if a single data file is available, an analysis can be made again.



The method contained in the data file is a copy of the method file that was used to acquire and analyze data. Therefore, when any method parameter in the data file opened via <Data Analysis> is modified, the method contained in the data file is modified rather than the method file.


 [Admin Manual]: [“4.1 Important File Concepts for Operation”](#)

1.2 Starting the LCMSsolution

This document assumes the following system configuration as an example to describe the procedure for an analysis:
High-pressure Gradient LCMS plus PDA (= Photo Diode Array) Detectors System

Pump	LC-10ADvp = 2 units
Auto sampler	SIL-10ADvp
Column oven	CTO-10A(C)vp
PDA detector	SPD-M10Avp
Mass spectrometer	LCMS-2010A


1 Check that the LC and MS units are On.



2 Check that nitrogen gas is sent to the MS unit.

3 Turn On the personal computer and peripheral devices to start Windows.

4 Enter your user ID to log on.





5 Double-click the [LCMSsolution] icon displayed on the Windows desktop. <LCMSsolution Launcher> will be started.



1.2 Starting the LCMSsolution

6 Select [Operation] menu.

7 Click the [Analysis] icon .
The <Login> screen will appear.




8 Select “Admin” and click the [OK] button.
The LCMS analysis program will be started with the <LCMS Analysis> main window displayed.

 [Admin Manual]: “2.4 Registering (Changing/Deleting) Users”, “2.5.2 Changing Passwords”





9 Click the [Data Acquisition] icon .



10 Check that “Ready” is displayed.

If “Not Connected” is displayed, properly complete <System Configuration>.

 [Operation Manual]: “14.5 Configuring System”



■ Description of <Data Acquisition> window

- Toolbar


Among the functions available on the Menu bar, the frequently used ones and the functions to directly control the analyzer are assigned to this bar.

- Assistant Bar

The icons to operate the application in accordance with the general analysis flow are assigned to this bar.

- Instrument Parameters

A pane is displayed showing the parameters for the system set up on <System Configuration>. Set those parameters for data acquisition.




2


Qualitative Processing (Single-run Analysis)

Set the parameters for the LC and MS units on the <Data Acquisition> window and then make an analysis. This document assumes an example of analysis under the following analytical conditions to specifically describe the procedure for the analysis.

Column	Shim-pack VP-ODS 150mm x 2.0mm i.d. 5µm (Equivalent to Shimadzu P/N 228-34937-94)
Mobile phase	Binary Gradient mode Pump A = Water, Pump B = Acetonitrile
Sample	Papaverine 0.5, 1, 5, 25, 50 ng/µL (Shimadzu P/N 225-06613-05)

2.1 Creating a new method file


- 1 Click the [New] button .
A new method file will be opened.



2.2 Setting the LC parameters

2.2.1 Detecting the auto sampler rack

- 1 Click [Advanced] button.
- 2 Select the [Autosampler] tab.
- 3 Click [Detect Rack] button.



2.2.2 Setting the LC parameters

[Operation Manual]: “4.2.1 Setting the LC Parameters”

- 1 Click [Normal] button
- 2 Select the [Simple Settings] tab.
- 3 Enter “6” min in [LC Stop Time].


If you click [Apply to all acquisition time] button after entering [LC Stop Time], [End Time] of all the detectors become the same.

- 4 Enter values for the pump parameters.

Mode	Binary gradient
T.Flow	0.2mL/min
B.Conc	45%

- 5 Enter “6” min in [End Time] of PDA.
- 6 Enter “40” °C for the oven temperature.

Be sure to enter a value in [Stop/End Time] (measurement end time) in steps 3 and 5.




2.2 Setting the LC parameters

■ Entering the gradient mode conditions

This document describes the procedure for setting up the pumps by assuming that liquid is sent in the gradient mode at a constant mixture ratio of the mobile phase.

To change the gradient mode conditions, perform the following steps:

- 1 Select the [LC Time Prog.] tab.
- 2 Enter values in [Time], [Module], [Action], and [Value] for the time program as shown on the right side.
- 3 Click [Draw curve] button.
The entered time program will be displayed as a graph.




■ Setting the pressure limit of a pump

If the column or the like is in an improper state, an error may occur because of exceeding pump's upper pressure limit. In this case, change the upper pressure limit by performing the following steps:

- 1 Click [Advanced] button.
- 2 Select [Pump] tab.
- 3 Enter "15" MPa in [P.Max].



The default value for [P.Max] is 10 MPa.



2.3 Setting the MS parameters



To set the MS (mass spectrometer) parameters, perform the following steps:

 [Operation Manual]: “4.2.2 Setting the MS Parameters”


- 1 Select the [MS] tab.
- 2 Enter “6” min in [Acquisition Time].
- 3 Select an event.
- 4 Set the parameters for the selected event.

Detector voltage	1.2kV
Measurement start m/z	200
Measurement end m/z	400

- 5 Click [DownLoad] button.
The instrument parameters will be transferred to the unit.
The dialog box will be opened allowing you to save the settings (method).
- 6 Enter “Method1.lcm” for the file name and click [Save].
The method file will be saved and the set parameters will be transferred to the unit.




Segment and Event





The LCMS-2010A provides the capability to allow you to change the analysis conditions in each specified time range during an analysis. The analysis conditions (a set of analysis conditions) in the specified time range are called a “Segment”. Multiple MS conditions may be specified for each segment and each of those conditions is called an “Event”. Additions of segments and events allow you to specify more complicated MS analysis conditions. This document assumes that an analysis is made under a single MS condition.

If multiple events are specified within the same segment, an analysis will be made under the condition specified for the event time and then the next event will occur. When the final event specified in the segment is finished, the first event will be resumed again. Thus, the cycle (Event#1 → Event#2 → Event#3 → Event#1... for Segment#1 in the above example) will be repeated for the time specified for the segment.

After the time specified for the segment has elapsed, similar operations will be performed for the next specified segment.

 If the “Polarity” (“Positive” or “Negative”) is changed, 400 msec is required for this change. This means that the time of the event after the polarity has been changed becomes shorter practically by 400 msec. Therefore, increase or decrease the event time as necessary.

 To add/delete any segment/event, right-click the appropriate segment/event in the event tree and select the desired option from the pop-up menu displayed.




2.4 Starting the operation of the instrument



Before starting an analysis, click the “Instrument Control bar” button at the top of the screen to start the operation of the analyzer. It will take about 20 minutes until the operation becomes stable enough.

2.4.1 Starting the control of the MS unit

- 1 Click the following five buttons: [Open/Close Nebulizing Gas], [CDL On/Off], [Heat Block On/Off], [IG On/Off] (= Ion Gauge On/Off), and [MS Detector On/Off].

The MS unit will start operating.

When an analysis is made using the APCI for the interface, the [APCI On/Off] button  is added.




- 2 For the LCMS-2010A, turn clockwise the knob for the drying gas controller to set the pressure.

For the LCMS-2010A-ESI: 0.1 MPa

For the LCMS-2010A-APCI: 0.02 MPa

Turn the knob clockwise.



2.4 Starting the operation of the instrument

2.4.2 Starting the operation of the LC unit

- 1 Click [Instrument On/Off] button.
The LC unit will start operating under the conditions specified in the method file.


HPLC Instrument




Item	Value	Units
Nebulizing Gas Flow	4.2	L/min
Nebulizing Gas Flow Monitor	45.1	L/min
CDL Temperature	250	C
CDL Temp. monitor	320	C
Heal Block Temperature	300	C
Heal Block Temp. monitor	300	C
Detector Voltage	1.20	kV
EG Vacuum	0.560	Pa
Flow	0.000	µL/min
Flow		
Pump Pressure	8.5	MPa
Oven Temperature	25.0	C
Maximum Temperature	65.0	C
Wavelength CDL	254	nm
Sample energy CDL	1200	nJ
Reference energy CDL	1544	nJ
CDL Max (Maximum)		
Detector Voltage		

2.4.3 Selecting a graph to be displayed in the <Chromatogram> view

The <Chromatogram> view allows you to specify the types and ranges of axes for the graph to be displayed.

 [Operation Manual]: “11.2 Customizing Windows”

- 1 Right-click anywhere on the graph and select the [Display Settings] menu.




- 2 Select the [MS] tab.
Enter values for m/z and other parameters for the mass chromatogram to be displayed.


- 3 Tick the check boxes on the 1st and 2nd rows.

- 4 Enter 340.15 on the 2nd row of the m/z column.

In this example, the mass chromatogram will be displayed according to TIC and $m/z = 340.15$.


- 5 Click [OK] button.


 To leave the <Display Settings> window open, click [Apply] button.




2.5 Acquiring data through a single-run analysis

To make a single-run analysis under the conditions specified in “2.2 Setting the LC parameters” and “2.3 Setting the MS parameters”, perform the following steps:

- 1 Click the [Single Start] icon . The <Single Run> window will be displayed.


 [Operation Manual]: “4.3 Starting a Single-run Analysis”



- 2 Enter “Sample1.lcd” for the data file name to be created.


- 3 Enter vial number “3” and injection amount “1”.

In this example, previously fill 5 ng/μL of papaverine into vial No. 3 of the auto sampler, and inject 1 μL from that vial.



- 4 Click [OK] button.


The single-run analysis will be started. After the [Acquisition Time] specified in the method file has elapsed, the analysis is finished automatically.





2.6 Performing qualitative processing on <MS Data Analysis>

2.6.1 Starting the <MS Data Analysis>


After the single-run analysis has been finished, perform data analysis as follows:

- 1 Click the [Data Analysis] icon . <MS Data Analysis> will be started. The last acquired data will be loaded and then displayed.

 [Operation Manual]: “5.1 Operation in the <MS Data Analysis> Window”




- 2 When the data file is first opened, only TIC is displayed in the <Chromatogram> View.



- 3 Dragging the cursor on each graph will allow you to enlarge that area. Right-clicking anywhere on each graph will allow you to select the [Initialize Zoom] or [Undo Zoom] option.

- 4 Clicking the [+] or [-] button will allow you to increase or decrease the level of the intensity axis.


- 5 Dragging the cursor on the splitter (frame) will allow you to change the aspect ratio of each view.



2.6.2 Displaying a mass spectrum


1 Double-click anywhere on the chromatogram.


The cut-out cursor will be moved to that time.
The mass spectrum for the cut-out cursor position in the <Chromatogram> View will be displayed in the <Spectrum> View.



Averaging the mass spectrum


Averaging the mass spectrum will allow you to obtain a clearer spectrum.

1 Click the [Average Spectrum] button  on the Toolbar.




2 Drag the cursor on the chromatogram to define the area you want to average.

The averaged spectrum in the defined time range (between 4.517 and 4.983 min in this example) will be displayed.



■ Performing subtractive processing of a mass spectrum


If the background mass spectrum is subtracted from the averaged spectrum, an even clearer spectrum can be obtained.

1 With the averaged spectrum displayed, click the [Average & Subtract Spectrum] button  on the Toolbar.

2 Drag the cursor on the chromatogram to define the area you want to subtract.

The spectrum obtained by subtracting the background will be displayed.

The information displayed above the spectrum graph indicates that the averaged spectrum for retention time between 3.400 and 4.227 min has been subtracted from that for retention time between 4.517 and 4.983 min.




■ Registering the averaged/subtracted spectrum in the “Spectrum Process Table”


If you register the averaged/subtracted spectrum in the spectrum processing table, you will be able to reproduce that spectrum easily on a later day.

1 Right-click anywhere on the spectrum graph and select [Register to Spectrum Process Table].

The averaged/subtracted mass spectrum will be registered.




Alternatively, it can also be registered by clicking the  button on the Toolbar.





2.6.3 Displaying a mass chromatogram

- 1 Double-click a mass spectrum peak.
A mass chromatogram will be additionally displayed in the <Chromatogram> View. The settings for the mass chromatogram are registered in the <Fragment Table> window.




Opening the <Fragment Table> window

- 1 Click the [Fragment Table] icon .
The <Fragment Table> window will be displayed.




Deleting the erroneously registered chromatogram

- 1 Remove a tick mark from the check box in the [Disp.] column on <Fragment Table> window.
- 2 Click [OK] button.
The window will be closed and the chromatogram will be hidden.



2.7 Performing peak integration (peak detection)

In this example, change the integration conditions in a single-run analysis and then perform peak integration again as follows:

- 1 Click the [Qualitative Peak Integration] icon .

The <Qualitative Peak Integration> window will be displayed.

- 2 Select the [Integration] tab.

- 3 Select "Detail" for the integration method.

If you select Auto (Area) or Auto (Height), peaks in the number close to the entered maximum number of peaks will be detected.

- 4 Enter "10" sec in Width.

If you specify the minimum width of peaks to be detected, the noise peak will be eliminated. Peaks will be detected to the extent that the half-width value is one forth the Width value.

- 5 Enter "1000" /min for the Slope value.

This is the parameter that determines the start and end points of the peak. When the absolute value of the gradient of the chromatogram becomes this value, the start and end points of the peak are determined there.

- 6 Click [OK] button.

The postrun will be carried out using the qualitative integration parameters you have set.

- 7 Click the [Qualitative Table] icon .


The <Qualitative Table> window will be displayed.

- 8 Select the [TIC] tab.

The integration result will be displayed.



The [Spectrum Process] tab allows you to check the registered averaged spectrum.



Peak	Ret. Time	Peak Start	Peak End	m/z	Area	Area%
1	4.759	4.433	5.467	TIC	104584013	1

2.7 Performing peak integration (peak detection)

Simple procedure for setting the integration parameters

Temporarily enter a little smaller values for Width and Slope and then double them, and see how peaks are detected*. In the example given in this document, first enter Width 10 and Slope 1000 and then Width 20 and Slope 2000.

* If the Width value is excessively increased, no minute noises will be detected as peaks.

If the Slope value is excessively increased, no moderate changes in the baseline will be detected as peaks.

Repeat the above steps and when the unnecessary peaks become undetectable, adopt the integration parameter at that point.


Checking data with <Data Explorer>


The LCMSsolution manages the data files, method files, batch files, and other related files in “Project Folders”.

<Data Explorer> allows you to manage the project of the LCMSsolution more effectively.

Project folders may be freely created, copied, or handled with <Data Explorer> of the LCMSsolution and the standard Explorer of Windows.

[Operation Manual]: “13.2 Managing Files Effectively”
[Admin Manual]: “6.1.1 Customizing Data Explorer Display Data”


1 Click the [Data Explorer] button .
This will toggle between displaying and hiding <Data Explorer>.



2 Change the display for each file type.

Double-clicking the file or dragging and dropping it to the window will allow you to load the file.

3 Right-click anywhere on the file icon.
A popup menu will appear.



Data Preview
The highlighted data file can be previewed.
Part of the sample information can also be checked.


Show File Info.
When “Detail” for [File View] is selected, the sample name and other additional information will also be displayed as the file information.


2.8 Printing out the analysis result

To print out the result of qualitative processing, perform the following steps.

2.8.1 Printing out a “Graph Image”

Print out the chromatogram and MS spectrum displayed on the screen as follows:



- 1 Click the [Print] button .
[Print Image] will be carried out.



■ Example of printing out a graph image


==== Shimadzu LCMSsolution Data Report ====

<Chromatogram>



<Spectrum>


Retention Time: 4.767(Scan#: 287)
Max Peak: 106 Base Peak: 340.10(4457918)
Spectrum: Averaged 4.750-4.783(286-288)
Background: Calc Polarity: Pos Segment: 1 - Event: 1




2.8.2 Selecting a layout for printing

<Data Report> allows you to print out a report image in the report format edited in the layout edit pane. In this example, load the preinstalled report format file “Sample1.lcr” to print out a graph image.


 [Operation Manual]: “10.2 Reprinting Data Processing Results”


- 1 Click the [Data Report] icon .
The data report will be displayed.



- 2 Select the [Report Format] tab with <Data Explorer>.

- 3 Drag and drop the file icon to the layout edit pane located on the right side.
The “Sample1.lcr” report format will be displayed.

- 4 Click the [Print] icon .
The report in the layout edit pane will be printed out.




2.8 Printing out the analysis result

Example of using the report format file for printing

Sample Information

```

Acquired by      : Admin
Date Acquired   : 6/19/2003 12:25:05 PM
Sample Type     : Unknown
Level#         : 0
Sample Name     :
Sample ID      :
ISTD Amount    : [1]=1 [2]=1 [3]=1 [4]=1 [5]=1
Sample Amount   : 1
Dilution Factor : 1
Vial#          : 2
Injection Volume : 1
Data File      : Sample1.lcd
Method File    : Method1.lcm
Original Method : C:\LabSolutions\Data\Project1\Method1.lcm
Report Format   :
Tuning File    : C:\LabSolutions\Data\Project1\ESI-1_5.lct
Processed by   : Admin
Modified Date  : 7/9/2003 10:43:18 AM
    
```



MS Peak Table TIC

Peak#	Ret. Time	Peak Start	Peak End	Area	Area%	Height	Height%	A/H	Mark	Name	ID#	Event
1	4.759	4.433	5.467	104584013	100.00	5761363	100.00	18.15		Papaverine	1	1-1
Total				104584013	100.00	5761363	100.00					

3

Quantitative Processing (Batch Analysis)

3.1 Creating a “Compound Table”


In the quantitative processing, the concentration of the compound contained in an “Unknown Sample” is calculated by creating a “Calibration Curve” with a “Standard Sample” of a known concentration, which contains the same compound as that being quantitatively analyzed.


In this example, inject 1 μL of a standard sample containing 0.5, 1, and 5 $\text{ng}/\mu\text{L}$ of papaverine to create a calibration curve. Simulate the quantitative processing to analyze 0.75 $\text{ng}/\mu\text{L}$ of papaverine as an unknown sample.

 [Operation Manual]: “5.5.2 Editing a “Compound Table””, “5.5.4 Using <Compound Table Wizard>”

3.1.1 Setting the quantitative parameters in <MS Data Analysis>

Set the quantitative parameters in the following steps using the papaverine data (Sample1.lcd) that has been loaded to <MS Data Analysis> in the previous chapter.

- 1 Click the  [Quantitative Parameters] icon.




- 2 Select the [Integration] tab.

- 3 Enter “100” /min for Slope.


In principle, enter a value equivalent to 1/2000 the targeted peak height.

If no peak is detected, reduce the Slope value by half.




3.1 Creating a “Compound Table”

- 4 Select the [Quantitative] tab.
- 5 Select “External Standard” for [Quantitative Method].
- 6 Enter “3” for [# of Calib. Levels].
- 7 Click [OK] button.





3.1.2 Creating a “Compound Table”


To complete the quantitative settings for each compound, set “Compound Table” to [Edit Mode].


- 1 Click [Edit] button  in the <Compound Table> View.
- 2 Enter values in the “Compound Table”.

Name	Type	m/z	Ret. Time	Conc. 1	Conc. 2	Conc. 3
Papaverine	Target	340.15	4.800	0.5	1	5


 If you click a peak in the <Chromatogram> View with the [Ret. Time] cell highlighted, the retention time for that chromatogram peak will be entered automatically.

 If you click a peak in the <Spectrum> view with the [m/z] cell highlighted, the m/z value for that spectrum peak will be entered automatically.

- 3 Click [View] button  .
The edited settings will be established.



■ Checking and saving the quantitative parameters/compound table


1 Click the [Peak Integration] icon .

2 Check for the identification mark (▼) on the chromatogram peak.

The identification mark is given to the identified peak.

The peak has the (↑) and (↓) marks at the starting and end points, respectively.



If the peak integration fails, adjust the Slope value among the integration parameters.

3 Check that the peak has been identified properly, and then click the [Apply to Method] icon .

The Save dialog box will be opened.

4 Check that “Method1.lcm” is selected for the file name, and then click [Save] button.

The method file will be overwritten.




3.2 Creating a SIM Table

The SIM (Selected Ion Monitoring) mode is the analysis mode that selects ion before data acquisition, and acquires the selected ions only. Therefore, the sensitivity is higher than the SCAN mode that acquires broader range of m/z values. In this example, use the mass number specified in “3.1 Creating a “Compound Table”” to change the parameters at data acquisition, so that the quantitative analysis can be made in the SIM mode at higher detection sensitivity.

1 Return to the <LCMS Analysis> - <Data Acquisition> screen.


2 Click the [MS] tab of the <Instrument Parameters> View.





3 Select “SIM” for the analysis mode.

The setting for the measured m/z value will be changed from entering a range to entering an individual m/z value.

4 For Ch1, enter “340.15” for the m/z value for papaverine in the compound table.





5 Click the [Save] button  to save the method file.




3.3 Creating a “Batch Table”

To make a batch analysis (continuous analysis), use the created method file to set up the batch table.

- Click the [Batch Processing] icon . The <Batch Table> window will be displayed. Create the batch table assigning the 1st to 3rd rows to standard samples and the fourth row to an unknown sample.




- Choose the [Edit]-[Table Easy Settings] menu. The <Table Easy Settings> window will be displayed.



- Select “New” for [Batch Table].
- Specify “Standard” samples.
Vial# : “1” - “3”
Data Filename : “Std01”

- Specify a “Unknown” sample.
Vial# : “4” - “4”
Data Filename : “Unk01”



- Click [OK] button. The 4-row batch table will be created.
- Tick the check box in the [Report Output] column and enter a file name in the [Report Format File] column.


Set only the fourth row for the unknown sample. In this example, specify the preinstalled report format file “Report1.lcr”.

	Vial#	Inj. Vol	Sample Name	Sample ID	Sample Type	Level#
1	1	1			1:Standard(1)	1
2	2	1			1:Standard	2
3	3	1			1:Standard	3
4	4	1			0:Unknown	0

Analysis Type	Method File	Data File	Report	Report Format File
IT QT MIT M	Method1.lcm	Std01.lcd	<input type="checkbox"/>	
IT QT MIT M	Method1.lcm	Std02.lcd	<input type="checkbox"/>	
IT QT MIT M	Method1.lcm	Std03.lcd	<input type="checkbox"/>	
IT QT MIT M	Method1.lcm	Unk01.lcd	<input checked="" type="checkbox"/>	Report1.lcr

7

3.3 Creating a “Batch Table”

 If the full path is not specified for any file name, the data will be created in the specified project folder.

The default values are given to the following items of the batch table. No modification is required of those values so far as the operations in this document are concerned.


- **Sample Type**

Clicking this column will display the <Sample Type> window shown on the right side.

Select a sample type from this window.

Select “Standard” for a sample to create/update a calibration curve or “Unknown” for a sample under quantitative analysis.

For the first standard sample to create a calibration curve, enable “Initialize Calibration Curve”.




- **Analysis Type**

Specify whether analytical processing is performed or not.

Clicking this column will display the <Analysis Type> window shown on the right side.

Tick the desired items.

For example, MIT (= Integration) shows that peak integration will be carried out and MQT (= Quantitative) indicates that quantitative calculation will be performed.



- **Level#**


Enter the level of a standard sample.


- **Report Output**


Ticking the check box will allow you to automatically print out the analysis result report.

- **Report Format File**


Clicking this column will display the <Select Report Format File> window shown on the right side. The analysis result report will be printed out in the report format specified here.








 [Operation Manual]: “9.3 Batch Processing Parameters”

8 Click the [Save] button  on the Toolbar.

9 Enter “Batch1.lcb” for the file name.



Entering data in the table cells


Type	Example	Description
Window popup type (for complicated settings)	 Click Here	When you click the button displayed to the right of the cell you have selected, the appropriate window pops up for you to enter data in that cell.
Drop-down list type (for selection from a list)	 Click Here	When you click the button displayed to the right of the cell you have selected, the available options are displayed in a drop-down list. Select the desired option from that list by clicking it.
Spin input type (for input of a specific value)	 Click here for increments Click here for decrements	When you click the upper or lower rectangle mark button displayed to the right of the cell you have selected, the stepped value assigned to that cell is increased or decreased. To enter any value other than the stepped values, directly enter it in the cell.
Check box type (for On/Off input)	 Click Here	Click the check box displayed on the cell to give or remove the tick mark.
Double-click type (for opening the file)	 Double-click a blank space	The data file or method file on the selected row of the batch table can be opened from the menu. Alternatively, the same operation can be performed by double-clicking a blank space in the cell.



3.4 Making a batch analysis


Using the batch table created in “3.3 Creating a “Batch Table”” make a batch analysis as follows.



1 Place the sample onto the autosampler.


Vial 1	Solution of 500 ppb papaverine (standard sample)
Vial 2	Solution of 1 ppm papaverine (standard sample)
Vial 3	Solution of 5 ppm papaverine (standard sample)
Vial 4	Unknown sample (to be determined) * In this example, a solution of 0.75 ppm papaverine is used as an unknown sample.

2 Click the [Batch Start] icon . During the batch analysis, <Batch Table> and the <Data Acquisition> window are simultaneously displayed in divided screens.

 To stop the batch analysis, click the [Stop] icon .

 If the batch analysis is paused, you may change the subsequent batch tables while continuing the analysis of the ongoing measurements.

 A snapshot can be performed to check the currently acquired data. To make a snapshot, click the [Snapshot] icon  on the [Acquisition] Assistant Bar during the analysis.





3 After the unknown sample has been analyzed, a report is output.

An example of printing out a report after batch analysis

Acquired by : Admin
 Date Acquired : 6/19/2003 12:01:05 PM
 Sample Type : Unknown
 Level# : 0
 Sample Name :
 Sample ID :
 ISTD Amount : [1]=1 [2]=1 [3]=1 [4]=1 [5]=1
 Sample Amount : 1
 Dilution Factor : 1
 Vial# : 10
 Injection Volume : 1
 Data File : Unk01.lcd
 Method File : Method1.lcm
 Original Method : C:\LabSolutions\Data\Project1\Method1.lcm
 Report Format : Report21cr
 Tuning File : C:\LabSolutions\Data\Project1\ESI-1_5.lct
 Processed by : Admin
 Modified Date : 8/8/2003 4:09:20 PM

Sample Information


MS Chromatogram
 C:\LabSolutions\Data\Project1\Unk01.lcd



ID# 1 m/z: 340.10
 Type: Target
 Name: papaverine
 Ret. Time: 4.743
 Area: 14816115
 Conc.: 0.707ppm

MS Calibration Curve

ID# 1 m/z: 340.10
 Name: papaverine
 Function: $f(x)=16,977,562.0X + 2,804,692$
 $Rr1=0.9999711$ $Rr2=0.9999423$
 Mean RF: 20,058,615.0 RF SD: 2,867,764 RF %RSD: 14.3




#	Conc.(ppm)	Mean Area
1	0.500	11591411.00
2	1.000	19447075.00
3	5.000	87729745.00

■ Printing out a summary of multiple results from batch analysis

After the batch analysis, print out a “Summary Report” (a simple report of more than one analysis result) as follows.

1 On <Report>, create a format for the summary report containing report items [MS Summary].

 [Operation Manual]: “10.3 Creating Report Files”



There are the following two types of summary report items:

- [Concentration]: The results of concentration, area, and height are displayed in a summary.
- [Compound]: The peak information such as concentration and column performance is displayed for each compound.

2 Enter [Summary Type] in <Batch Table>.

Specify “Summary Start” for the top of the data to be output to the summary report, “Summary Run” for the data to be included in the summary report, and “Summary End” for the data on the final line to be included in the summary report.

3 Enter [Summary Report Format File].


Enter a file name to the right of the cell in which you have specified “Summary Start”.



For example, if you complete the following settings, the summary report including the data “Tutorial_Unk01.lcd”, “Tutorial_Unk02.lcd”, and “Tutorial_Unk04.lcd” will be printed out in the format “Summary1.lcr” when the batch analysis is finished.

	Analysis Type	Method File	Data File	Summary Type	Summary Report Format
1	IT QT MIT MQT	Tutorial_Method.lcm	Tutorial_Unk01.lcd	Summary Start	Summary1.lcr
2	IT QT MIT MQT	Tutorial_Method.lcm	Tutorial_Unk02.lcd	Summary Run	
3	IT QT MIT MQT	Tutorial_Method.lcm	Tutorial_Unk03.lcd	None	
4	IT QT MIT MQT	Tutorial_Method.lcm	Tutorial_Unk04.lcd	Summary End	

②
③




4 Click the [Batch Start] icon .

The batch analysis will be made.

After the batch analysis has been finished, the specified summary report file is printed out.

■ An example of printing out a summary report

MS Summary(Compound)
Tutorial Unk01.lcd



[MS] ID1 Compound Name:papaverine

Title	Sample Name	Sample ID	Retention Time	Area	Height	Concentration
Tutorial Unk01.lcd			4.743	11077384	835152	0.746
Tutorial Unk02.lcd			4.737	10886185	846829	0.735
Tutorial Unk04.lcd			4.718	10418105	807638	0.706
Average			4.733	10793891	829873	0.729
%RSD			0.276	3.142	2.425	2.840
Maximum			4.743	11077384	846829	0.746
Minimum			4.718	10418105	807638	0.706
Std. Dev.			0.013	339191	20122	0.021


4 Data Analysis

4.1 Checking a “Calibration Curve”


To check and modify the “Calibration Curve” that has been created using the data on the standard sample analyzed in Chapter 3, use <MS Calibration Curve>.


1 Select the [Method] tab of <Data Explorer> displayed in <LCMS Postrun>.

2 Double-click the method file “Method1.lcm”.




3 Select [Close & Open] button in the selection dialog box.


 Since the method file “Method1.lcm” is loaded by <Data Acquisition> in Chapter 3, temporarily close the file.





4 <MS Calibration Curve> will be displayed.

5 Check to see whether all of three area values are registered.

-  If the area value is 0, no peak integration has been performed. Adjust the Slope value and then carry out peak integration again.



-  To change the “Slope” value: Click the [Quantitative Parameters] icon  on the [MSCalibration] Assistant Bar and then change the “Slope” value on the [Integration] tab.

-  To perform “Peak Integration” again: Click the [Peak Integration for All Data] icon  on the [MSCalibration] Assistant Bar to carry out peak integration.


6 Click the [Save] button on the Toolbar.

The modified method file will be saved.

4.2 Checking the quantitative calculation result of an unknown sample

Using the <MS Data Analysis> window, check the data analysis result of the unknown sample analyzed in Chapter 3 as follows.


1 Click the [Data] tab of <Data Explorer> displayed on <LCMS Postrun>.



2 Double-click the data file “Unk01.lcd” that has been obtained by analyzing the unknown sample.



<MS Data Analysis> will be displayed with the data file “Unk01.lcd” loaded.


3 If the calibration curve has been changed in “4.1 Checking a “Calibration Curve””, then drag and drop the data file “Method1.lcm” from <Data Explorer> - [Method] tab to <MS Data Analysis>. The “Calibration Curve Information” will be imported.



4 Check that the identification mark (▼) is displayed on the chromatogram peak. If the mark is not displayed, the peak integration has not been completed successfully. Adjust the Slope value and then carry out the peak integration again.

 To change the “Slope” value: Click the [Quantitative Parameters] icon  on the [MS Data] Assistant Bar and then change the “Slope” value on the [Integration] tab.

 To perform “Peak Integration” again: Click the [Peak Integration] icon  on the [MSData] Assistant Bar to carry out peak integration.


5 Click the [Save] button  on the Toolbar.


The reanalyzed data files will be saved.


4.3 Loading a batch file to the “Quant Browser”


Using of the “Quant Browser” (= Quantitation Browser) allows you to easily re-analyze multiple data.


4.3.1 Displaying the quantitative result from the batch file

- 1 Click the [LCMS Browser] icon  on <LCMSsolution Launcher>. <LCMS Browser> will be displayed.


 [Operation Manual]: “8.1 Browsing the Quantitative Calc. Results at a Time”




- 2 Click the [MS Quant Browser] icon  on the [LCMS Browser] Assistant Bar.




- 3 Drag and drop the file icon “Tutorial_Batch.lcb” from the [Batch] tab of <Data Explorer> to <Quant Browser>. All of the sample data (“Tutorial_Std01.lcb” through “Tutorial_Std03.lcb” and “Tutorial_Unk01.lcb” through “Tutorial_Unk05.lcb”) will be loaded.


 Alternatively, the data may be loaded by selecting multiple data file from <Data Explorer> and then simultaneously dragging and dropping them.




4.3 Loading a batch file to the “Quant Browser”

4 Click the compound table.
The quantitative result of the compound on specified row will be displayed.

 To delete the data file, right-click the <Quantitative Result> View and then choose [Delete] from the menu displayed.

 The calibration curve for the above compound will also be displayed.

5 Check the chromatogram in the <Chromatogram> View.



4.3.2 Setting the integration parameters again to retry peak integration


The sample data consists of the quantitative data obtained using three-points absolute calibrations.


However, it shows that the data processing of the standard sample in a low concentration range is failed due to improper integration parameters.


In this example, the data on the highlighted 1st row of the <Quantitative Result> View indicates a failure when the file has been loaded.


If you check the area value, you will find that it is zero.

If you also check the <Chromatogram> View, you will see that the peaks have not been integrated.





1 Click the  [Quantitative Parameters] icon.


 [Operation Manual]: “8.2 Making a Postrun Analysis of Multiple Data”



- 2 Click the [Integration] tab.
- 3 Enter “100” /min for the [Slope] value.
If the value is too large, enter a smaller value.
- 4 Click [OK] button.



- 5 Click the [Peak Integration for All Data] icon  to retry the peak integration.
The peaks will be detected.
The three-point calibration curve will be displayed.



The proper quantitative value has been obtained.

Data Filename	Area	Height	Conc. (ppm)	Std. Conc.
Tutorial_Std01.lcd	0.00000
Tutorial_Std02.lcd	15,230,436	1,152.3	1.00000	1.00000
Tutorial_Std03.lcd	80,746,261	4,833.3	5.00000	5.00000
Tutorial_Unk01.lcd	11,077,384	835.2	0.74644
Tutorial_Unk02.lcd	10,986,185	846.3	0.73477
Tutorial_Unk03.lcd	10,791,516	841.4	0.72899
Tutorial_Unk04.lcd	10,419,105	807.4	0.70619
Tutorial_Unk05.lcd	10,493,710	799.0	0.71080
Average	21,377,657	1,445.3	1.37531
%RSD	122,715,247	103,747.5	116.458396
Maximum	80,746,261	4,833.3	5.00000
Minimum	10,419,105	799.0	0.70619
Std. Dev.	26,233,644.2351	1,499,355.0	1.601667



Data Filename	Area	Height	Conc. (ppm)	Std. Conc.
Tutorial_Std01.lcd	11,591,411	737.6	0.51755	0.50000
Tutorial_Std02.lcd	19,447,075	1,253.8	0.98026	1.00000
Tutorial_Std03.lcd	87,729,745	4,931.3	5.00219	5.00000
Tutorial_Unk01.lcd	14,816,115	927.6	0.70749
Tutorial_Unk02.lcd	14,840,655	956.6	0.70893
Tutorial_Unk03.lcd	14,803,802	952.0	0.70676
Tutorial_Unk04.lcd	14,238,404	913.9	0.67346
Tutorial_Unk05.lcd	14,004,370	897.3	0.66439
Average	23,943,947	1,446.3	1.24513
%RSD	108,018,348	97,859.2	122.349813
Maximum	87,729,745	4,931.3	5.00219
Minimum	11,591,411	737.6	0.51755
Std. Dev.	25,863,856.0442	1,415,338.38	1.523413

4.3 Loading a batch file to the “Quant Browser”

■ Files handled by the Quant Browser

The <Quant Browser> is an application that reanalyzes multiple data using the same method file for data processing. Files are loaded in accordance with the following rules:

- **Method file**

Load a method file from the [Method] tab of <Data Explorer>.

If you do not specify a method file, the method file of the first loaded data will be loaded automatically.

If the loaded method file contains calibration curve information, the data file for the standard sample used to create that calibration curve will be loaded.

- **Data file**


Load a data file or data files from the [Data] tab of <Data Explorer>.


The use of the Toolbar buttons allows you to display the data for each sample type.



4.4 Printing out a summary report from the Quant Browser

<Quant Browser> has the “Summary Report” capability to report all of the loaded data as follows.

- 1 Click the [Summary Report] icon . The image for each compound in the table will be printed out.




■ An example of printing from the quantitation browser

=== Shimadzu LCMSSolution Quant. Browser Report ===


[MS] ID1 Compound Name:papaverine						
Title	Sample Name	Sample ID	Ret Time	Area	Height	Conc.
Tutorial_Std01.lcd			4.682	11591411	737667	0.518
Tutorial_Std02.lcd			4.717	19447075	1253846	0.980
Tutorial_Std03.lcd			4.736	87729745	4931305	5.002
Tutorial_Unk01.lcd			4.743	14816115	927676	0.707
Tutorial_Unk02.lcd			4.737	14840655	956688	0.709
Tutorial_Unk03.lcd			4.699	14803802	952027	0.707
Tutorial_Unk04.lcd			4.718	14238404	913809	0.673
Tutorial_Unk05.lcd			4.708	14084370	897395	0.664
Average			4.717	23943947	1446302	1.245
%RSD			0.444	108.018	97.859	122.350
Maximum			4.743	87729745	4931305	5.002
Minimum			4.682	11591411	737667	0.518
Std. Dev.			0.021	25863856	1415339	1.523


4.5 Using the Data Browser


Using of the data browser allows you to display multiple data files in various types of formats as follows.

- 1 Click the [Browser] icon  from <LCMSSolution Launcher>. <LCMS Browser> will be started.

 [Operation Manual]: "8.3 Listing Multiple Data"





- 2 Click the [Data Browser] icon .






- 3 Open a data file (multiple data files may be selected) by dragging and dropping it.


A window will pop up for the user to select the data type displayed, whether to replace data or add a cell, and the direction of adding that cell in the latter case.



 A maximum of 64 cells (8 x 8) may be displayed.

 If you click the focus pin  located in the upper right corner of each cell so that the pin is displayed as  in multiple cells, then the display of those cells will be changed interlocking with each other.

For example, if both cells for an MS chromatogram and a PDA spectrum are “pinned”, then double-clicking the MS chromatogram will display the PDA spectrum for that time.




5

Exiting the LCMSsolution


5.1 Existing the LCMSsolution

1 Click **X** in the upper right corner of the screen.


Alternatively, you may select the [Exit] menu located at the bottom of the [File] menu to exit the LCMSsolution programs.



2 When using <LCMS Analysis>, click [OK] button on the confirming dialog box.



3 For <LCMS Analysis>, the <Shutdown> window is displayed. Give the tick mark to all the check boxes.




4 Click [OK] button. The shutdown procedure will be started.

5 Click [No] button in the confirmation dialog box.

For any file that has not been saved, the confirming dialog box is displayed to prompt you to confirm whether the file must be saved when you exit the LCMSsolution.

6 All the LCMSsolution programs will be terminated with Windows shut down.



Index

Guide: A
 Operation: B
 Administration: C

Numerics

21 CFR Part 11 C1, C51

<3D Image> B66

A

Absolute Retention Time Method B59

Acquisition Mode B14

Acquisition Type B14

Acrobat B148, C49, C57

Distiller C50, C57, C59

File Attachment Tool C52

PDFWriter C57

Action (Batch Table) B97, C18

Add/Specify PC C56

[Add Calibration Level] button B93

<Add Group> C11

[Add/Remove Programs] icon B143

<Add/Remove Programs Properties> B143

[Add Rows with Selected Data File] button B95

<Add Rows with Selected Data Files...> B95, B124

<Add User> C12

Admin B5, C3, C5

Adobe

Acrobat B148, C49, C57

File Attachment Tool C52

Acrobat Reader B136, B148, C50

Installation Procedure B148

Adobe PostScript Driver C58

Agent C2

<Agent Registration Settings> C54

AIA(ANDi)

Data Management C51

Output (Batch Table) C17

All-in-One Structure C25, C51

[All users must use password] check box C9

Analysis Mode A28

<Analysis Type> A30

Analysis Type (Batch Table) A30, B97

[Analyze] icon B57, B59, B61, C31

APCI A11

[Apply to all acquisition time] button B19

[Apply to Method] icon A27, B53, B56, B62, B81, C32

Archiver C51

Arithmetic Operations B82

ASCII Output (Batch Table) C17

Assistant Bar

[Real Time]

[Tuning]icon B7

Assistant Bar A5, B5, B113, B114

[Acquisition]

[Batch Processing] icon A29

[Data Analysis] icon A15, A32, B12

[Single Start] icon B16

[Snapshot] icon A32, B12

[Stop] icon B20

[Acquisition]

[Single Start] icon A14

[Batch]

[Batch Start] icon A32, A34, B94, B96, C38

[Pause/Restart] icon B94, B100

[Setting] icon C17, C18

[Settings] icon C33

[Stop] icon A32, B20, B100

[Wizard] icon B92

[Data Report]

[Preview] icon B103

[Print] icon A23, B103

[LC Data]

[Analyze] icon B57, B59, B61, C31

[Apply to Method] icon B53, B56, B62, B81, C32

[Data Report] icon B53, B56, B103

[LCMS Brower]

[Data Browser] icon B87

[LC Quant Brower] icon B83

[MS Quant Brower] icon B83

[LCMS Brower]

[Data Browser] icon A44

[MS Quant Brower] icon A39

[LC Quant Brower]

[Data Analysis Parameters] icon B85

[Modify Calibration Curve] icon B85

[Summary Report] icon B86

[MSCalibration]

[Peak Integration for All Data] icon A37

- [Quantitative Parameters] icon A37
- [MS Data]
 - [Apply to Method] icon A27, B53, B56, B62, B81, C32
 - [Data Report] icon A23, B53, B56
 - [Fragment Table] icon A18, B29
 - [Peak Integration] icon A27, A38, B49
 - [Qualitative Peak Integration] icon A19
 - [Qualitative Table] icon B28, B31, B33
 - [Quantitative Parameters] icon A25, A38, B38, B40
 - [Wizard] icon B42
- [MS Data Analysis]
 - [Qualitative Table] icon A19
- [MS Quant Brower]
 - [Peak Integration for All Data] icon A41
 - [Quantitative Parameters] icon A40
- [MS Quant Browser]
 - [Modify Calibration Curve] icon B51, B85
 - [Peak Integration for All Data] icon B86
 - [Quantitative Parameters] icon B85
 - [Summary Report] icon A43, B86
- [PDA Data]
 - [Analyze] icon B57, B59, B61, C31
 - [Apply to Method] icon B53, B56, B62, B81, C32
 - [Data Report] icon B53, B56, B103
 - [Spectrum Index] icon B66
 - [Wizard] icon B42
- [Postrun]
 - [Batch Processing] icon B95, C38
 - [LC Data Analysis] icon B57, B59, B61, C31
 - [Library Editor] icon B36
 - [MS Calibration Curve] icon B51
 - [Report Format] icon B105, B106
- [Real Time]
 - [Batch Processing] icon B91
 - [Data Acquisition] icon A4, B22
 - [System Configuration] icon B15, B17, B22, B149, C21
 - [System Control] icon B24
 - [Tuning] icon B8, B9
- [Report]
 - [Preview] icon B105, B109
 - [Print] icon B104
- [Tuning]
 - [Auto Tuning Condition] icon B9
 - [Manual Tuning View] icon B9
 - [Start Auto Tuning] icon B7
 - [Tuning Result View] icon B8
- [Assistant Bar] button B12
- <Audit Trail> C28, C30
- Audit Trail
 - Data File C22
 - Method File C22
 - Setting C21, C22, C25, C26
- Audit Trail Log
 - Merge C10
 - Print B111

- Auto Purge (Batch Table) B97
- Auto Tuning B7, B10
- [Auto Tuning Condition] icon B9
- Average (Bracket Calibration) B98
- [Average Spectrum] button A16
- Average Spectrum B32
- [Average & Subtract Spectrum] button A17

B

- Background Data File (Batch Table) B99
- Baseline Check (Batch Table) B97
- [Baseline Check] button B11
- Baseline Chromatogram B99
- Base Period B21, B22
- [Base Shift] menu B64
- Batch Analysis A25, A29, A34, B48
- Batch File A2, B6, B126, B127, B129, B130, C23
 - Prohibit overwriting C10
- Batch Processing A34
 - Action B97, C18
 - AIA(ANDi) Output C17
 - Analysis Type A30, B97
 - ASCII Output C17
 - Auto Purge B97
 - Background Data File B99
 - Baseline Check B97
 - Custom Parameters B97
 - External Program Run C18
 - Level# A30, B51, B63
 - Parameter B97
 - Progress Control C18
 - QA/QC Criteria C17, C32
 - QA/QC Function C33
 - Report Format File A30, B102
 - Report Output A30
 - Run Mode B97, B99
 - Sample Type A30, B96
 - Start A32
 - Stop A32
 - Summary Report Format File A34, B97
 - Summary Type A34, B97
 - System Check B97, C17
 - Tuning File B97
 - User Program B97, C18
- [Batch Processing] icon A29, B91, B95, C38

- <Batch Queue> B100
 - Batch Schedule B96
 - [Batch Start] icon A32, A34, B94, B96, C38
 - <Batch Table> A29, B113, B129
 - <Add Rows with Selected Data Files...> B95, B124
 - Create A29
 - <Fill Down...> B93, B123
 - Print B111
 - <Settings>
 - [Bracket] tab B98
 - [Data Filename] tab B91, B96, C44
 - [QA/QC] tab C33
 - [Shutdown] tab B98
 - [Startup] tab B98
 - <Table Easy Settings> A29, B123
 - <Batch Table Wizard> B92, B102, B123
 - Blank Peak Rejection B99
 - Browsing File B84
- ## C
-
- <Calibration Curve> B84, B113, B129
 - Calibration Curve A36, B43
 - Create B48
 - <Calibration Curve> View B65
 - Calibration Point
 - Add B48, B51, B63
 - Delete B51, B63
 - [Cancel Edit] menu B39, B62
 - CDL Heater A11
 - Cell B117
 - Multiple Select B121
 - Title Cell B117
 - Value Entry A31
 - Cell Data
 - Clear B121
 - Copy B121
 - Cut B121
 - Paste B121
 - Value Enter B122
 - Cell (Data Browser) A45
 - Adjust Layout B90
 - Display Link B88
 - Display Settings B89
 - <Change Analysis Time> B20
 - [Change Database] button C55
 - <Check Raw Data> C29
 - <Check the Program Files> C37
 - <Choose Destination Location> B139
 - Chromatogram
 - Export B81
 - Extract B69, B70
 - Operation B82
 - Registration B71
 - Chromatogram File
 - CLASS-LC10 B133
 - <Chromatogram> View A13, A40, B27, B65, B68, B69, B70, B75
 - <Display Settings> A13, B11
 - Magnify B64
 - <Properties> B11
 - Shift B64
 - Stack/Overlapping mode B71, B75
 - Stretch/Contraction B64
 - Undo B64
 - CLASS-8000
 - Data File B133
 - Method File B133
 - CLASS-Agent C2, C51, C54, C57
 - CLASS-LC10
 - Chromatogram File B133
 - Data File B133
 - Method File B133
 - CLASS-VP
 - Data File B133
 - Method File B133
 - Clipboard B55
 - Column Performance Test C31
 - Compound Finder B47
 - “Compound Table” A25, A40, B39, B44, B51, B62, B63, B83
 - [Conc.] column B41
 - [ISTD Group] column B41
 - [Name] B62
 - [Ret. Time] B62
 - [Type] B62
 - [Type] column B41
 - <Compound Table> View B25, B39, B62, B65
 - [Edit] button A26
 - [Param’s] tab B39, B44, B47, B62, B67, B80, B84
 - [Results] tab B49, B84
 - [View] button A26
 - <Compound Table Wizard> B42

<Confirm File Deletion> B144
Consumable C35
 Reset B151
<Contour> View B65, B68, B69, B70
[Control Panel] B143
<Convert LCsolution Data file to ASCII file> B132
[Create filename automatically with...] check box B96
<Create New Project (Folder)> B128
<Create SIM Table> B54
Criteria B100
Custom Parameters (Batch Table) B97

D

DAO B137
DAO Installation B145
Data Access Objects B137
<Data Acquisition> A5, A28, B11, B113, B129
[Data Acquisition] icon A4, B22
<Data Analysis> B25, B113, B129
[Data Analysis] icon A15, A32, B12
<Data Analysis Parameters>
 [Identification] tab B59
 [Integration] tab B21, B57, B68
 [Library] tab B76
 [Multi Chrom] tab B67
 [Performance] tab C31
 [Purity] tab B74
 [Quantitative] tab B61
 [UV Spectrum] tab B78
 Width B21
[Data Analysis Parameters] button B57, B59, B61, B68, B74, B76, B78, C31
[Data Analysis Parameters] icon B85
<Database Maintenance> C55
 [Change Database] button C55
 [Network] check box C55, C56
<Data Browser> A44, B26, B87, B129
[Data Browser] icon A44, B87
<Data Explorer> A20, A23, A39, B5, B113, B127
 Display Mode C44
 Display Settings C43
 Docking View C44
 Extension tab C48

Floating View C44
[Data Explorer] button A20, C44
<Data Explorer Properties> C48
Data File A2, A42, B6, B56, B104, B126, B127, B129, C23
 All-in-One Structure C25, C51
 Attachment to PDF C52
 Audit Trail C22
 Automatic Filename Creation B91, C44
 CLASS-8000 B133
 CLASS-LC10 B133
 CLASS-VP B133
 LCMSsolution Ver.2.x B133
 Long-Term Storage C1, C49
 Method in Data File C28
 Prohibit overwriting C10
 Rollback C28
[Data File Properties...] menu C26
<Data Report> A23
Data Report C25, C60
 Print B103
[Data Report] icon A23, B53, B56, B103
Data Rollback C28
<Data Tree> View B51, B52, B63
Detect Rack A7, B13
Deviation Rejection of Baseline B99
Differential Chromatogram B99
[Display All] button B65
[Display Contour View] button B65
[Display Multi Chromatogram Table...] menu B67
[Display Purity View] button B65
<Display Settings> A13, B26, B89
Distiller C50, C57
Double-Click B114, B129
[DownLoad] button A9
[Download] button B16
Drying Gas Controller A11

E

Easy Login mode C10
<Edit Compound Information> B37
Electronic Record C1, C48
Electronic Signature C1

E-Mail function
 Disable C10
Event A10, B14
Event Log C13
<Exit Setup> B140
Export
 Chromatogram B81
 Method File B56, B81, C25
 Report Format File B53, B56
 Spectrum B81
Exposure time B155
Extension B6, B127
External Program Run (Batch Table) C18

F

FDA 21 CFR Part 11 C1, C51
File
 Convert B132
 Copy B130
 Create C19
 Delete B131
 Double-Click B114, B129
 Drag B129
 Drop B129
 Extension B6, B127
 History Information B126, B131
 Move B130
 Open B125, B129
 Double-Click B114, B129
 Drag & Drop B129
 Rename B131
 Save B126
<File New> C19
<File Properties> C22, C26
 [Audit Trail] tab C22
 [Sample Info] tab C26
 [Used Files] tab C26
<File Search> C45
 Comment C45
 Creation Date C45
 Filename C45
 Sample Information C45
 Search Result C46
 Type C45
 Update Date C45
File Type B6

<Fill Down...> B93, B123
Focus Pin button A45, B87
<Fragment Table> A18, B33
Fragment Table B29
[Fragment Table] icon A18, B29
Function Check C38
 Data Installation C40
 Evaluation Criteria C39
 Inspection Report C42

G

GLP/GMP C1
Good Laboratory Practice C1
Good Manufacturing Practice C1
Grouping B39, B43, B62
<Group Property> C11
GxP Support Function C2

H

Hardware Validation C34
<Header/Footer> B109
Heat Block A11
Help B135
Hit List B35, B36
'H/W Administrator' Rights C4

I

Identification Parameters B59
"Identification Results Table" B49
Import
 Method File B56, C25
 Report Format File B53, B56
Infusion Analysis B17
[Infusion Control] bar B17
 Operation B18
[Initialize Calibration Curve] button B93
Installation Procedure B137
Install CD-ROM B137
<Install Complete> B140
Install Disc B137

- <Instrument> B149
 - [Instrument Control] bar B113
 - [Baseline Check] button B11
 - [Shutdown] button B23
 - [Startup] button B23
 - [System Check] button C34
 - [Instrument Monitor] bar B12, B113
 - [Instrument Monitor] button B12
 - [Instrument On/Off] button A12
 - <Instrument Parameters> C27
 - Instrument Parameters A5, C24
 - LC A7, B13, B19
 - MS A9
 - Instrument Parameter Settings B11
 - <Instrument Parameters> View A7
 - [Auto Purge] tab B13
 - [Autosampler] tab A7
 - [Data Acquisition] tab B21, B22
 - [Sampling] list box B21, B22
 - Display Mode B13
 - [DownLoad] button A9
 - [Download] button B16
 - [LC Time Prog.] tab B13, A8
 - [MS] tab A28, B14, B19
 - [Normal/Advanced] button A7
 - [Normal]/[Advanced] button B13
 - [Pump] tab A8
 - [Simple Settings] tab A7, B13, B19
 - Integration Parameters A40, B57
 - <Integration Time Program> B31, B58
 - Internal Standard Method B40
 - Internal Standard Substance B41
 - Ion Gauge A11
- ## L
-
- <LabSolutions LCMSsolution Setup> B137, B142, B145
 - [Standard Install] button B138
 - LabSolutions PDF Printer C60
 - Layout
 - Save B26
 - Layout File B90
 - Save B90
 - lcb A2, B6, B127
 - <LC Calibration Curve> B63
 - lcd A2, B6, B127
 - <LC Data Analysis> B55, C47
 - [LC Data Analysis] icon B57, B59, B61, C31
 - <LC Data Comparison> B64, B87, B113, B129
 - LC Instrument Parameters B13, B19
 - lcm A2, B6, B127
 - <LC Maintenance Utility> C35
 - LC Maintenance Utility B151
 - <LCMS Analysis> A4, A28, B7, B113
 - <LCMS Analysis Editor> B5
 - [LCMS Analysis Editor] icon A1, B3
 - [LCMS Analysis] icon A1, A4, B3, B7, B149
 - <LCMS Browser> A39, A44, B83, B87, B113
 - [LCMS Browser] icon A1, A39, A44, B3, B83, B87
 - <LCMS Postrun> B113
 - [LCMS Postrun] icon A1, B3
 - <LCMSsolution Launcher> A3
 - [LCMS Analysis] icon A1
 - <LCMSsolution Lancer>
 - [LCMS Analysis] icon A4
 - <LCMSsolution Launcher> B3, B7, B91, B95, C8, C9, C11, C12, C13, C14, C38
 - [LCMS Analysis Editor] icon A1, B3
 - [LCMS Analysis] icon B3
 - [LCMS Browser] icon A1, A39, A44, B3, B83, B87
 - [LCMS Postrun] icon A1, B3
 - [Log Browser] icon B4, C13
 - [System Policy] icon B4, C9
 - [User Administration] icon B4, C11, C12, C14
 - <LCMSsolution Launchre>
 - [LCMS Analysis] icon B149
 - LCMSsolution Ver.2.x
 - Data File B133
 - Method File B133
 - <LCMS Tool Launcher> B151, B153
 - lcq B84
 - [LC Quant Browser] icon B83
 - lcr A2, B6, B103, B127
 - LC Stop Time B19
 - Level# (Batch Table) A30, B51, B63
 - <Library Editor> B113, B129, B130
 - [Library Editor] icon B36

Library File B126, B129, B130, C24
 Edit B37
Library Search B34, B77, B87, B90
[Library Search] menu B35
Library Search Result
 Print B111
<Library Search Results> B34, B35
Library Type B79
<License Installer> B141, B142
License Key Disk B141
License Uninstallation B142
[Lock] menu C16
<Log Browser> C13
 Print B112
[Log Browser] icon B4, C13
<Login> A4, B5, B149, C10, C55, C56
Login
 Attempt Limit C7, C9
 Lockout Duration C7, C9
Login/Logout Condition Log C30
[Login with Windows Username] check box C7, C9
<Log of Data Acquisition> C27
LSSMntDB.exe C55
LSSMntPC.exe C56
Luminous Energy B156
lyt B90

M

<Maintenance of Instrument Information> C56
Manual Identification B50, B52
Manual Peak Identification B85, B86
Manual Peak Integration B31, B50, B52, B85, B86
Manual Tuning B9
[Manual Tuning View] icon B9
Margin B109
Mass Chromatogram A18
Mass Number Calibration B9
Mass Pattern Matching B46
Mass Spectrum A16
Mass-to-Charge Ratio B9, B27

Maximum Number of Hits B34
Maximum Password Age C7, C9
Max Plot B67
Menu bar B5
Message Log C30
Method
 Print B111
'Method Developer' Rights C4
Method File A2, A42, B6, B51, B63, B83, B126, B127, B129, C24
 Audit Trail C22
 CLASS-8000 B133
 CLASS-LC10 B133
 CLASS-VP B133
 Create A6
 Export B53, B56, B62, B81, C25
 Import B53, B56, C25
 LCMSsolution Ver.2.x B133
 Prohibit overwriting C10
 Save A9, B62
[Method File Properties...] menu C22
Method in Data File C28
Method Validation C31
<MIC Table> B29
Minimum Similarity Level B34
Modify Calibration Curve B86
[Modify Calibration Curve] icon B51, B85
MS
 Start and Stop B24
 System Control B24
<MS Calibration Curve> A36, B51, B129
[MS Calibration Curve] icon B51
<MS Data Analysis> A15, A38, B25
<MS Fragment Table> B89
MS Instrument Parameters B14, B19
<MS Library Editor> B36
[MS Quant Browser] icon A39, B83
MS Stop Time B19
<MS Tuning> B129
[Multi-Chrom Table] button B67, B68, B70
[Multi User] check box C7, C9
Multi-User Mode C6
m/z B9, B27

N

Nebulizing Gas A11
Negative A10
[Network] check box C55, C56
[New] button A6, B36, B79, B105
Noise Spectrum B74
[Normal/Advanced] button A7

O

Online Help B13, B18, B135
Online Manual B136
'Operator' Rights C4
<Options>
 [Auto Copy] tab C51
 [PDF Export] tab C49
<Output Window> B5, B113, C30
 [LogFile] tab C30
 [Message Log] tab C30
Overlap (Bracket Calibration) B98

P

<Page Setup> B109
Page Size B109
Part 11 C1, C51
Password B5, C5, C12, C14, C15
 Change C14
 Duplication C9
 Force C9
 Maximum Age C7, C9
 Mean Decryption Time C7
 Minimum Length C7, C9, C12, C14
Pattern Matching B46
[Pause/Restart] icon B94, B100
PC restart B144
<PDA Data Analysis> B65
<PDA Utility> B153, C35
PDA Utility B153
 <Instrument Information> B157
 <Wavelength Calibration> B155
 <Wavelength Check> B154
PDF File C51

 File Information C50
PDF Format C49
 Output Setting C57
PDFWriter C50, C57
Peak Detection A19
[Peak Integration by ID] button B86
[Peak Integration for All Data] icon A37, A41, B86
[Peak Integration for All IDs] menu B86
[Peak Integration] icon A27, A38, B49
Peak Integration Parameters A19
Polarity A10
Policy C7, C9
 Level Setting C9
Portable Document Format C49, C57
Positive A10
Postrun Batch Processing B48, B52, B96
 Start B96
Post Search B36
Pre-filter B79
<Preview> B109
Preview C43, C46
[Preview] button B105, B109
[Preview] icon B103, B105, B109
Print
 Audit Trail Log B111
 Batch Table B111
 Data Report B103
 Graph Image A21
 Image B111
 Image Extraction from PDF C53
 Library Search Result B111
 Log B112
 Method B111
 Postrun Batch Processing B102
 Program Alteration Check B112
 Raw Data Integrity Check B112
 Real Time Batch Processing B102
 Single Run B101
 Spectrum Table B111
 System Check B112
 System Configuration B111
[Print] button A21, B105
[Print] icon A23, B103, B104
Print Spectrum Index B111

Profile B14

Program Alteration Check C37
 Print B112

Progress Control (Batch Table) C18

Project B127, C51
 Create B128
 Delete B128

Project Folder A20

<Properties> B26
 [Sample Info.] tab B49

<Properties> (<Chromatogram> View) B11

<Properties> (File) C22, C26
 [Audit Trail] tab C22
 [Sample Info] tab C26
 [Used Files] tab C26

<Properties> (<System Configuration>) B21
 [Base Period] B21

<Purity> View B65
 [Peak Profile] tab B75
 [Peak Purity] tab B75

Q

QA/QC Function B100

<QA/QC>
 [Criteria] tab C32
 [General] tab C32

QA/QC Criteria (Batch Table) C17, C32

QA/QC Function (Batch Table) C33

[QA/QC Parameters] menu C31, C32

Qualitative Analysis B33

[Qualitative Layout] button B65

<Qualitative Parameters>
 [Library] tab B34
 [Spectrum Process] tab B30

<Qualitative Peak Integration> A19

[Qualitative Peak Integration] icon A19

<Qualitative Table> B28
 [Peak Table] tab B33
 [Spectrum Process] tab B28, B34

[Qualitative Table] icon A19, B28, B31, B33

<Quant Browser> A39, B26, B51, B83, B113, B129

<Quantitation> View B49, B50

Quantitative Calculation Method B43

Quantitative Calculation Result Check B49

[Quantitative Layout] button B65

<Quantitative Parameters> B38, B85
 [Identification] tab B45, B46
 [Integration] tab B38
 [Quantitative] tab B40

Quantitative Parameters B61

[Quantitative Parameters] icon A25, A37, A38, A40, B38, B40

[Quantitative Paramters] icon B85

<Quantitative Result> View A40, B83

R

Raw Data
 Integrity Check C29
 Print B112

Real Time Batch Processing B48, B94
 Additional Analysis B100
 Pause B94, B100
 Restart B94
 Resume B100
 Start B94

[Real Time Instrument Control] bar A11, B113
 [APCI On/Off] button A11
 [CDL On/Off] button A11
 [Heat Block On/Off] button A11
 [IG On/Off] button A11
 [Instrument On/Off] button A12
 [MS Detector On/Off] button A11
 [Open/Close Nebulizing Gas] button A11
 [Zeros Detector] button B11

<Reference Ion Ratio> B45

Reference Ion Ratio B52

Reference Ion Ratio Identification B45

[Register Current Spectrum as Standard] menu B44

Registered Spectrum B77

[Register Spectrum to Library] menu B35

[Register Standard Spectrum] menu B80

[Register to Multi-Chrom Table] button B71

[Register to Spectrum Process Table] button A17

[Register to Spectrum Table] button B72

Registration
 Data to Agent C54
 Rights Group C11
 User C5, C12

- Reinstallation B142
 - Relative Retention Time Method B41, B59, B60
 - <Report> A34, B113, B129, B130
 - Report
 - Multiple number of data files B104
 - Report Format A23, A34, B92
 - Template B103
 - Report Format File A2, A23, B6, B103, B104, B105, B126, B127, B129, B130, C23
 - Export B53, B56
 - Import B53, B56
 - Prohibit overwriting C10
 - Report Format File (Batch Table) B102, A30
 - [Report Format] icon B105, B106
 - Report Format in Data File B56, C25
 - Report Item A34, B110
 - Format B106, B107
 - Layout B108
 - Paste B106
 - Resize B107
 - [Report Item] bar B113
 - Report Output (Batch Table) A30
 - Report Type (QA/QC) C32
 - <Reset Lamp Time> C35
 - Right
 - 'System Administration' C4, C11, C12
 - Use Analytical Instrument C12
 - Rights Group C4, C7, C12
 - 'H/W Administrator' C4
 - 'Method Developer' C4
 - Name C11
 - 'Operator' C4
 - Registration C11
 - [Rollback to Original Data] menu C28
 - Row# Column B117
 - Run Mode (Batch Table) B97, B99
- ## S
-
- Sample Information B55, C43
 - Sample Injection
 - Repetition B93
 - Sample Loading Method B17
 - Sample Schedule B94
 - <Sample Type> A30
 - [Add Calibration Level] button B93
 - [Initialize Calibration Curve] button B93
 - Sample Type (Batch Table) A30, B96
 - Sample Type (QA/QC) C32
 - [Sampling] list box B21, B22
 - Sampling Rate B21, B22
 - [Save Batch in Data File As...] menu C28
 - [Save] button A31, A28, A37, A38, B93, B96
 - [Save Library File As...] menu B79
 - [Save Method As...] menu C25
 - [Save Report Format File As...] menu B103
 - Scan B14
 - Security Level C7
 - Security Support Function
 - Current User Display C15
 - Login Authentication C15
 - Restrict the Operations Right C16
 - Window Lock C16
 - Segment A10, B14
 - <Select Baseline> B31, B50
 - <Select Batch Table> C28
 - <Select Components> B146
 - <Select Data Type> A44, B87
 - Sequence (Bracket Calibration) B98
 - <Set Compound Information> B45
 - [Setting] icon C17, C18
 - <Settings>
 - [Bracket] tab B98
 - [Data Filename] tab B91, B96, C44
 - [Create filename automatically with...] check box B96
 - [QA/QC] tab C33
 - [Shutdown] tab B98
 - [Startup] tab B98
 - Shift
 - Left-Right B64
 - Upper-Lower B64
 - <Show Standard Spectrum> B46
 - [Show Window] menu B114
 - <Shutdown> A46, B23
 - Shutdown B98
 - [Shutdown] button B23
 - SIM A28, B14, B29, B54

- Similarity B77
- Similarity Identification B80
- SIM Table A28
 - Create B54
- <Single Run> A14, B101
- Single-run Analysis B48, A6, A14
- [Single Start] button B16
- [Single Start] icon A14, B16
- Single-User Mode C6
- Snapshot
 - Update B53, B81
- [Snapshot] icon A32, B12
- Software Validation C37
- Spectrum
 - Background Processing B30
 - Export B81
 - Extract B69, B70
 - Operation B82
 - Registration B72
- Spectrum-based Identification B44
- Spectrum File C24
- Spectrum Index
 - Print B111
- [Spectrum Index] icon B66
- Spectrum Process Table A17, B28, B33, B35
- <Spectrum Table> B73
- Spectrum Table
 - Print B111
- [Spectrum Table] button B72
- <Spectrum> View B27, B35, B65, B68, B69, B70, B77
 - Extracted/Registered mode B72
- [Standard Install] button B138
- Standard Sample B48
- [Standard Sample] button B9
- Standard Spectrum B46
- [Start Auto Tuning] icon B7
- <Startup> B23
- Startup B98
- [Startup] button B23
- Statistical Calculation B86
- Status Information C24
- [Stop] button B20
- [Stop] icon A32, B20, B100
- Stop Time
 - LC B19
 - MS B19
- Stretch/Contraction
 - Left-Right B64
 - Upper-Lower B64
- Subtract Spectrum B32
- Summary Report Format File (Batch Table) A34, B97
- [Summary Report] icon A43, B86
- Summary Report Output (Quant Browser) A43
- Summary Type (Batch Table) A34, B97
- System Administration C3
 - History C13
 - Log Browser C13
 - Network C56
 - Policy C7
- 'System Administration' Right C4, C11, C12
- "System Administration" Tools C8
- <System Check> C34
 - Print B112
- System Check B8
- System Check System Configuration
 - Criteria B150
- System Check (Batch Table) B97, C17
- [System Check] button C34
- <System Configuration> A4, B11, B150, C27, C30
 - <MS Detector> B17
 - Print B111
 - <Properties> B21
 - [Base Period] B21
- System Configuration
 - Audit Trail C21
 - Refer C27
- <System Configuration Audit Trail Settings> C21
- System Configuration File C21, C24
- [System Configuration] icon B15, B17, B22, B149, C21
- [System Control] icon B24
- <System Policy> C7, C10
 - [All users must use password] check box C9
 - [Login with Windows Username] check box C7, C9
 - [Multi User] check box C7, C9
- [System Policy] icon B4, C9
- System Suitability Test C31

T

[Tab Control] button B116

Table Column B117

- Display Condition Change B118
- Select B117
- Width Change B117

<Table Easy Settings> A29, B123

Table Row

- Add B119
- Delete B119
- Display Quantity Change B120
- Select B119

<Table Style> B26, B39, B40, B62, B80, B97, B118, B120

Template

- Create from Template C19
- Save as a Template C20

Title bar B5

Title Cell B117

Toolbar A5, B5, B113

- [Add Rows with Selected Data File] button B95
- [Assistant Bar] button B12
- [Average Spectrum] button A16
- [Average & Subtract Spectrum] button A17
- [Data Analysis Parameters] button B57, B59, B61, B68, B74, B76, B78, C31
- [Data Explorer] button A20, C44
- [Display All] button B65
- [Display Contour View] button B65
- [Display Purity View] button B65
- [Instrument Monitor] button B12
- [Multi-Chrom Table] button B67, B68, B70
- [New] button A6, B36, B79, B105
- [Peak Integration by ID] button B86
- [Preview] button B105, B109
- [Print] button A21, B105
- [Qualitative Layout] button B65
- [Quantitative Layout] button B65
- [Register to Multi-Chrom Table] button B71
- [Register to Spectrum Process Table] button A17
- [Register to Spectrum Table] button B72
- [Save] button A28, A31, A37, A38, B93, B96
- [Single Start] button B16
- [Spectrum Table] button B72
- [Standard Sample] button B9
- [Stop] button B20
- [Tab Control] button B116
- [Update Snapshot] button B53, B81

Total time of lamp illumination B156

<Tuning> B7

- [Manual Tuning] B9
- [Tuning Result] B8

Tuning B24

- Auto Tuning B7, B10
- Manual Tuning B9
- Standard Sample B9

Tuning File B129, C24

Tuning File (Batch Table) B97

[Tuning] icon B8, B9

[Tuning]icon B7

[Tuning Result View] icon B8

U

Uninstallation Procedure B142

Unknown Sample B48

[Update Snapshot] button B53, B81

<Update Wizard> B47

[Use Password] check box C12

User C3

- Registration C5, C12

<User Administration>

- [Groups] tab C7, C11
- [Password] tab C14
- [Users] tab C12

User Administration

- Description C3
- Group C3
- Instrument C3
- Password C3
- Policy C9
- User ID B5, C3, C12
- User Name C3, C12

[User Administration] icon B4, C11, C12, C14

<User Authentication Mode Select> B139

User ID B5, C3, C12

<User Information> B138

User Name C3, C12

User Program (Batch Table) B97, C18

<User Property> C12, C14

- [Use Password] check box C12

<UV Library Editor> B36, B79

UV Library File B76, B79, B80
<UV Library Search Results> B77
UV Spectrum B80

V

Vacuum Pump B24

Validation

Hardware C34
Method C31
Software C37

View

Magnified B25, B55, B65
Size Change B115
Standard B25, B55, B65

W

<Wavelength Settings> B67
Width (Data Analysis Parameters) B21
Window/Band Method B60
Window Display Format B115
(Docking/Floating view) B116
[Wizard] icon B42, B92

Z

[Zeros Detector] button B11

